




Fluoromount-G® Anti-Fade

Cat. No.	Format	Size
0100-35	Liquid mounting medium	20 mL



Human pancreatic carcinoma cell line MIA PaCa-2 was stained with Mouse Anti-Human CD44-BIOT (SB Cat. No. 9400-08; right) followed by Streptavidin-TXRD (SB Cat. No. 7100-07), DAPI, and mounted with Fluoromount-G® Anti-Fade (SB Cat. No. 0100-35).

Description

Fluoromount-G[®] Anti-Fade is an aqueous-based compound for mounting slides used in microscopy applications in which the preservation of fluorescence intensity is critical. This mounting medium provides a semi-permanent seal for storage of slide preparations and since Fluoromount-G[®] Anti-Fade is water-soluble, the coverslip may be removed by submerging the slide in a PBS solution until the coverslip is loosened.

Applications

IHC-FS IHC-PS ICC

Fluoromount-G® Anti-Fade has been shown to be compatible with the following reagents based on in-house analysis –

Alexa Fluor® 488 CY3 Rhodamine (TRITC)
Alexa Fluor® 555 CY3.5 Texas Red®

CY2 Fluorescein (FITC)

Suggested Protocol

- Prepare samples and perform staining according to established procedures
- Following final step, remove excess liquid around the sample and add 1 drop (~50 μL) of Fluoromount-G[®] Anti-Fade for an area of 22 x 22 mm
- Mount coverslip and remove excess mounting medium
- Analyze

Handling and Storage

• Fluoromount-G® Anti-Fade is supplied as a 20 mL solution. Avoid exposure to light and store at 2-8°C.

Fluoromount-G[®] is a registered trademark of Southern Biotechnology Associates, Inc. Alexa Fluor[®] and Texas Red[®] are registered trademarks of Thermo Fisher Scientific, Inc. and its subsidiaries Cy[®] is a registered trademark of GE Healthcare.

For Research Use Only. Not for Diagnostic or Therapeutic Use.

TB0100-35 24-Oct-19