



U.S. Corporate Headquarters
400 Valley Rd.
Warrington, PA 18976
1(800) 523-2575 / (215) 343-6484
1(800)343-3291 fax
info@polysciences.com

Polysciences Europe GmbH
Handelsstrasse 3
D-69214 Eppelheim, Germany
+(49) 6221-765767
+(49) 6221-764620 fax
info@polysciences.de

Polysciences Asia-Pacific, Inc.
2F-1, 207 Dunhua N. Rd.
Taipei, Taiwan 10595
(886) 2 8712 0600
(886) 2 8712 2677 fax
info@polysciences.tw

TECHNICAL DATA SHEET 1024

page 1 of 3

Transporter™ 5 Transfection Reagent

OVERVIEW:

Transporter™ 5 transfection reagent is the ready-to-use form of a proprietary linear polyethylenimine (PEI) derivative. Transporter™ 5 condenses DNA into positively charged complexes which readily enter the cell by endocytosis. The Transporter™ 5-DNA complex is highly resistant to endosomal degradation, allowing high transfection efficiency.

Transporter™ 5 transfection reagent effectively introduces DNA into HEK-293, CHO, COS, HeLa, insect (Sf9 and Sf21) and other eukaryotic cell lines. Transporter™ 5 is most compatible with plasmids up to 135,000 nucleotides in size, although larger plasmids can be used with some success.

REAGENTS INCLUDED:

- **Transporter™ 5 transfection reagent:** 1-5 x 1 mL vial of reagent containing 1 µg/µL proprietary linear PEI derivative. 1 mL is sufficient for transfecting a total of 250 µg of DNA, or about 60-120 transfections in 6 well plates.

REAGENTS TO PREPARE:

- **Diluent:** A 150 mM NaCl solution prepared with sterile, WFI/Cell Culture Water. A solution of 25 mL (containing 0.2912 g NaCl) is sufficient for 1 mL of Transporter™ 5 transfection reagent.

PROTOCOL FOR ADHERENT CELLS:

Preparation:

- Plate cells 18 to 24 hours before transfection.
- Use an appropriate number of cells in seeding solution to obtain a cell monolayer with 60-80% confluence. These qualities will provide the most optimal conditions for transfection. See Table 1 for suggested guidelines.
- Note: High serum levels inhibit the efficacy of Transporter 5™. In most cases, low serum levels ($\leq 5\%$) will produce the highest transfection efficiency.

Table 1- Guidelines for seeding adherent cell culture vessels.

Culture Vessel	Total Surface area (cm ²)	Cells in seeding solution
96-well plate	0.3	[1.2 – 2.4] x 10 ⁴
48-well plate	1.0	[4.0 – 8.0] x 10 ⁴
24-well plate	1.9	[0.8 – 1.6] x 10 ⁵
12-well plate	3.5	[1.5 – 3.0] x 10 ⁵
6-well plate	9.6	[4.0 – 8.0] x 10 ⁵
35 mm dish	9.6	[3.5 – 7.0] x 10 ⁵
60 mm dish	21	[0.9 – 1.8] x 10 ⁶
100 mm dish	58	[2.2 – 4.4] x 10 ⁶
T75 flask	75	[3.0 – 6.0] x 10 ⁶
T175 flask	175	[0.7 – 1.4] x 10 ⁷

Sample Transfection Protocol (Single well in 6-well plate):

1. 1 to 2 hours before transfection, exchange growth media in single well with 3 mL fresh growth media containing 2% serum.
2. Prepare Transporter™ 5-DNA transfection mixture (order is critical):
 - i. To 300 µL diluent in polypropylene tube add 2 µg plasmid DNA.
 - ii. Briefly mix/vortex solution.
 - iii. Add 8 µL Transporter™ 5 to mixture. (1:4 DNA/Transporter)
 - iv. Vortex for 5 seconds.
 - v. Let solution sit for 20 minutes in hooded environment to allow Transporter™ 5-DNA complexes to form.
 - vi. Mix solution gently by pipetting up and down 3 times.
3. Add Transporter™ 5-DNA transfection mixture to well.

Note: The preceding protocol can easily be scaled up or down by adjusting the volume of the Transporter™ 5-DNA transfection mixture in diluent to 10% of the overall culture volume. Ensure that the DNA/Transporter™ 5 ratio is 1:4. See Table 2 for recommended amounts of reagents for various culture vessels.

Incubation

- Following addition of Transporter™ 5-DNA transfection mixture to well, return wells to incubator.
- Typically, recombinant protein is detectable at 36-48 hours after transfection. Maximal expression is usually observed 72-96 hours after transfection.

Table 2 - Reagent quantities for transfection in a variety of culture vessels.

Culture Vessel	Culture Volume (mL)	Plasmid DNA (µg)	Diluent (mL)	Transporter™ 5 (µL)
6-well plate, single well	3	2 – 4	0.3	8 – 16
35 mm dish	3	2 – 4	0.3	8 – 16
60 mm dish	5	6 – 12	0.5	24 – 48
100 mm dish	10	12 – 24	1.0	48 – 96
T75 flask	15	18 – 36	1.5	72 – 144
250 mL shake flask	50	50 – 100	2.5	200 – 400

PROTOCOL FOR SUSPENSION CELLS:

Preparation:

- 2 to 3 hours before transfection, seed cells at 1.0×10^6 per mL of culture.

Sample Transfection Protocol (50 mL culture in 250 mL shake flask)

1. Prepare Transporter™ 5-DNA transfection mixture (order is critical):
 - i. To 2.5 mL of diluent in polypropylene tube, add 50 µg of plasmid DNA.
 - ii. Briefly mix/vortex solution.
 - iii. Add 200 µL of Transporter™ 5 to mixture.
 - iv. Vortex solution for 5 seconds.
 - v. Let solution sit for 20 minutes in hooded environment to allow Transporter™ 5-DNA complexes to form.
 - vi. Gently mix solution by pipetting up and down 3 times.

2. Add entire transfection solution to 25 mL of cell suspension culture.
3. Return cell suspension to incubator.

Incubation

- Shake for 2 to 3 hours in incubator, then add 25 mL of fresh culture medium. Return to incubator.
- Typically, recombinant protein is detectable at 36-48 hours after transfection. Maximal expression is usually observed 72-96 hours after transfection.

ORDERING INFORMATION

Catalog No.	Description	Size(s)
26008	Transporter™ 5 Transfection Reagent	(1 mL), 26008-5 (5 x 1 mL)

RELATED PRODUCTS

26253	Cell Culture Water	500 mL, 1000 mL, 6 x 500 mL, 6 x 1000 mL
25626	Cell Culture Bag, 50 mL with Filment	15 bags per box
25628	Cell Culture Bag, 250 mL with Filment	10 bags per box
25631	Cell Culture Bag, 1 L with Filment	5 bags per box
25634	Cell Culture Bag, 2 L with Filment	5 bags per box
25627	Cell Culture Bag, 50 mL with 2 Ports	15 bags per box
25630	Cell Culture Bag, 250 mL with 2 Ports	10 bags per box
25633	Cell Culture Bag, 1 L with 3 Ports	5 bags per box
25636	Cell Culture Bag, 2 L with 3 Ports	5 bags per box
25629	Cell Culture Bag, 250 mL with Filment & 2 Ports	10 bags per box
25632	Cell Culture Bag, 1 L with Filment & 3 Ports	5 bags per box
25635	Cell Culture Bag, 2 L with Filment & 3 Ports	5 bags per box

TO ORDER

In The U.S. Call: 1(800) 523-2575 • (215) 343-6484

In The U.S. Fax: 1(800) 343-3291 • (215) 343-0214

In Germany Call: +(49) 6221-765767

In Germany Fax: +(49) 6221-764620

In Asia Call: (886) 2 8712 0600

In Asia Fax: (886) 2 8712 2677

Order online anytime at www.polysciences.com